AVALON LABORATORIES

AVALON ELITE®

Multi-Port Venous Femoral Catheter

Instructions for Use
Multi-Port Venous Femoral Catheter

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





F

General Statement

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Note: Any medical techniques described in this document are meant for informational purposes only and should not replace the expertise of the physician and/or medical protocols.

Indication for Use

The Avalon Elite® Multi-Port Venous Femoral Catheter with Introducer is indicated to drain central venous blood via the femoral vein during extracorporeal procedures for up to six hours.

Device Description

The device is supplied with an introducer to facilitate placement into the vasculature by normal access techniques. Both the catheter and introducer are radiopaque and include depth marks. The catheter is wire-reinforced. The introducer is designed to follow a pre-positioned standard 0.038" (0.97mm) guide wire (not included).

The device is not MRI compatible.

How Supplied

The device is packaged in a peel pouch and carton. The device is STERILE (EtO), NON-PYROGENIC unless opened or damaged.

Instructions for Use:

Note: Consult Figures A-F while reviewing the following instructions.

1. Inspection: Check packaging and device for damage prior to use.

Warning: Do not use if package is open or damaged. Follow an aseptic technique to remove the catheter and introducer from the pouch.

2. Sizes: Selection should be based on patient anatomy and flow requirements.

Warning: Incorrect sizing can result in vascular damage and inadequate flow performance.

3. Device Assembly: Remove device from its pouch. Remove catheter and introducer from their protective sheaths. Straighten introducer and insert it into the catheter body. Fully engage the cap onto the connector. Adjust the introducer extension as needed.

4. Insertion and Placement: Always follow standard percutaneous or surgical techniques to place the device. Catheter depth marks are referenced from the tip of the catheter.

Warning: Incorrect insertion can cause damage to the vessels and/or heart structures.

Caution: Always use the introducer (included) over a 0.038" (0.97 mm) guide wire (not included).

Caution: Verify the advancement, positioning, and placement of the guide wire using appropriate imaging technology.

Caution: Verify the advancement, positioning, and placement of the catheter with pre-loaded introducer using appropriate imaging technology.

Caution: Never advance, position, or place the catheter without the introducer as this may kink and damage the catheter.

Caution: Blood leakage from the proximal drainage holes may occur during catheter insertion until these holes are fully advanced into the vessel.

Caution: Ensure adequate vessel size when the device is placed in order to allow for distal perfusion of the limb.

5. Securing the Device: Secure the device to the patient following standard medical practice.

Caution: A suture tied directly around the wire-reinforced section can cut, kink, or damage the catheter.

6. Introducer Removal and Clamping: After verifying correct placement, slowly and carefully remove the guide wire and then the introducer.

Caution: Do not clamp the catheter until the introducer tip has cleared the clamp tube section.

Caution: Blood loss may occur if the introducer is removed from the cap prior to clamping.

Caution: Do not clamp the wire-reinforced section as this will cause permanent damage to the catheter.

Caution: Do not clamp within 2 cm of the rigid barbed connector as this may cause permanent damage to the connector.

7. Device Removal: Clamp the non wire-reinforced tube section and then carefully remove the catheter from the vessel. Close the entry site following standard medical practice.

8. Disposal: Dispose of used products in accordance with established hospital protocol for biohazard waste.

Warning: The device is intended for single use only. Do not resterilize or reprocess as this can damage the device. Do not reuse as this could injure the patient and/or cause communication of infectious agents.

Contraindication

The physician is responsible for determining any patient conditions or perfusion complication which may prohibit the safe use of the device. The device should only be used as indicated.

Potential Complications

Possible complications are those normally associated with extracorporeal procedures, anticoagulation, and large size venous femoral catheterization.

Additional Warnings and Precautions

- Only physicians trained and experienced with femoral venous catheterization in conjunction with extracorporeal procedures should use this device.
- Do not use without adequate systemic anticoagulation.
- Standard medical practice for catheter insertion site care should be followed to reduce incidence of infection.
- Do not modify the device in any way.
- Do not use alcohol or alcohol-based fluids for lubrication or cleaning.

Avalon Laboratories, LLC
Limited Warranty and Disclaimer

Avalon Laboratories LLC (“Avalon”) warrants that each product (“Product”) it manufactures will be free from defects in materials or workmanship prior to its “Use By Date.” Avalon’s sole liability shall be limited to replacement of, or credit for, any Product returned to Avalon, upon determination by Avalon, in its sole discretion, to have been defective in materials or workmanship. To qualify for this limited warranty, the Product must be returned to Avalon within forty-five (45) days after usage of the Product and must not have been altered or subjected to misuse, abuse or accident. In addition, the Product must be used in accordance with the labeling and instructions for use provided with the Product.

Except as expressly set forth above, Avalon disclaims all other warranties, with respect to the Product, including but not limited to any IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Avalon shall not be liable to any person or entity for any medical expenses or damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. AVALON IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT.

This Agreement shall be governed by and construed in accordance with the laws of the State of Delaware without regard to any conflict of law or choice of law rules or provisions.

In no event shall the United Nations Convention on the International Sale of Goods apply to this Agreement or to any transactions under the Agreement.

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AVALON LABORATORIES, LLC