


## Bi-Caval Dual Lumen Catheter

### Catheter Dimensions



Description	A	B			C		D		E		F		G		H			
	Connector Size	Body Diameter			Insertable Length		Proximal Insertable Body to Infusion Port		SVC to Infusion Port Length		Infusion Port to Tip Length (Lead Length)		SVC to Tip Length		Lead Diameter			
Units	in	Fr.	mm	in	cm	in	cm	in	cm	in	cm	in	cm	in	Fr.	mm	in	
Description Size	Order Code																	
13 Fr.	10013	1/4	13	4.3	0.17	10	4.0	8	3.0	1.7	0.7	2.8	1.1	4.5	1.8	11	3.7	0.15
16 Fr.	10016	1/4	16	5.3	0.21	13	5.0	9	3.6	3.1	1.2	4.0	1.6	7.1	2.8	14	4.7	0.19
19 Fr.	10019	1/4	19	6.4	0.25	20	7.75	15	5.8	4.3	1.7	5.7	2.3	10.1	4.0	16	5.3	0.21
20 Fr.	10020	3/8	20	6.6	0.26	29	11.5	21	8.3	5.3	2.1	9.4	3.7	14.7	5.8	17	5.7	0.22
23 Fr.	10023	3/8	23	7.7	0.30	29	11.5	21	8.3	5.3	2.1	9.4	3.7	14.7	5.8	20	6.7	0.26
27 Fr.	10027	3/8	27	9.0	0.35	29	11.5	21	8.3	5.3	2.1	9.4	3.7	14.7	5.8	24	8.0	0.31
31 Fr.	10031	3/8	31	10.3	0.41	29	11.5	21	8.3	5.3	2.1	9.4	3.7	14.7	5.8	27	9.0	0.35

Tolerances may vary, specifications for reference use only and are subject to change.