




Multi-Port Venous Femoral Catheter

Catheter Dimensions



Description	A	B		C		D		E		F		G		H			
	Connector Size	Insertable Length		Proximal Insertable Body to Side Port		First Proximal Side Port to Second Port		Side Port to Side Port Length		Tip to Side Port Length		Tip to Proximal Port Length		Catheter Outer Diameter			
Units	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	Fr.	mm	in	
Description Size	Order Code																
20 Fr.	11020	3/8	61.0	24.0	29.5	11.6	7.6	3.0	3.8	1.5	4.8	1.9	31.5	12.4	20	6.7	0.263
22 Fr.	11022	3/8	61.0	24.0	29.5	11.6	7.6	3.0	3.8	1.5	4.8	1.9	31.5	12.4	22	7.3	0.288
24 Fr.	11024	3/8	61.0	24.0	29.5	11.6	7.6	3.0	3.8	1.5	4.8	1.9	31.5	12.4	24	8.0	0.315
26 Fr.	11026	3/8	61.0	24.0	29.5	11.6	7.6	3.0	3.8	1.5	4.8	1.9	31.5	12.4	26	8.7	0.341
28 Fr.	11028	3/8	61.0	24.0	29.5	11.6	7.6	3.0	3.8	1.5	4.8	1.9	31.5	12.4	28	9.3	0.367

Tolerances may vary, specifications for reference use only and are subject to change.