

# **AVALONELITE**® Vascular Access Kit

*Patent Pending*

## The **AVALONELITE** Advantage

- + Extended dilator tip is designed for smooth vessel insertion and guide wire tracking
- + Large range of dilator sizes accept standard 0.038" guide wires
- + Non-coated stainless steel guide wire displays 10cm incremental depth markings
- + 3mm J-tip guide wire end can be temporarily straightened to aid in placement

## Specifications

### 100cm Guide Wire Kit


- + 10Fr., 12Fr., 16Fr., 20Fr., 24Fr. Dilator
- + 0.038" 100cm Guide Wire with 10cm depth mark increments
- + 18Ga. Introducer Needle
- + #11 Blade Mini Scalpel
- + 10cc Syringe

### 210cm Guide Wire Kit

- + 10Fr., 12Fr., 16Fr., 20Fr., 24Fr., 26Fr., 30Fr. Dilator
- + 0.038" 210cm Guide Wire with 10cm depth mark increments
- + 18Ga. Introducer Needle
- + #11 Blade Mini Scalpel
- + 10cc Syringe

Learn more at:

[www.avalonlabs.com/vak](http://www.avalonlabs.com/vak)



Order Code	Guide Wire Length	Order Unit
12100	100 cm	5/Carton
12210	210 cm	5/Carton